

Fitting Guide for Raphael GP contact lenses.

Cylinder power of Cornea (D)	Diameter of the lens (mm)		
	Less than 8.90	9.00 – 9.50	Over 9.60
	First choice BC in „mm” – related to the flatter „K readings” (on „K”)		
0.00	0.05 mm steeper or on “K”	0.05 mm steeper	0.10 – 0.15 mm flatter
0.25 – 0.75	0.05 – 0.10 mm steeper	„On K”, or 0.05 mm steeper	0.05 – 0.10 mm flatter
0.75 – 1.75	0.10 – 0.15 mm steeper	0.05 – 0.10 mm steeper	„On K” or 0.10 mm flatter
1.75 – 2.75	0.15 – 0.20 mm steeper	0.15- 0.20 mm steeper	„On K” or 0.10 mm steeper
Over 2.75	0.25 – 0.30 mm steeper	0.20 – 0.25 mm steeper	0.10 – 0.20 mm steeper

The „on K” means the fatter K readings!

Of course always need a controll under „BURTON LAMP” and the eye need to be tinted by with „FLOURETS strip”.

Have to modify the fitting if needed after this controll.

Acuity Contact Lenses

A revolution in optical engineering